

3M Press-On™ Prism Availability

PRISMS

- Common Usage:
Vision Therapy, e.g.,
- Strabismus
 - Amblyopia
 - Diplopia



Sold one lens per package

| | |
|-------------------|--------------------|
| 1.00 ^Δ | 10.00 ^Δ |
| 2.00 ^Δ | 12.00 ^Δ |
| 3.00 ^Δ | 15.00 ^Δ |
| 4.00 ^Δ | 20.00 ^Δ |
| 5.00 ^Δ | 25.00 ^Δ |
| 6.00 ^Δ | 30.00 ^Δ |
| 7.00 ^Δ | 35.00 ^Δ |
| 8.00 ^Δ | 40.00 ^Δ |
| 9.00 ^Δ | |

PLUS SPHERES

- Common Usage:
- Low Vision
 - Recreational, e.g.,
divers masks, swim/ski goggles



Sold one lens per package

Aspheric Plus

| | |
|-------|--------|
| +0.50 | +4.00 |
| +1.00 | +6.00 |
| +1.50 | +8.00 |
| +2.00 | +10.00 |
| +2.50 | +12.00 |
| +3.00 | +14.00 |
| +3.50 | +16.00 |

GENERAL PRODUCT INFORMATION:

Press-On Optics:

- Apply smooth side to the back surface of one or both lenses of the patient's eyewear.
- May cover the entire lens or any region of a conventional lens.
- Can be easily cut with an ordinary pair of scissors or razor blade.
- Adhere to the lens by using plain water (hydrostatic pressure)—no glue, chemicals or residue to harm the carrier lenses.
- Add no noticeable weight—only 1mm thick, regardless of the prescription strength. Soft and pliable.
- Call Signet Armorlite Technical Services at 800 759-0075 if you need help with Press-On application.

LITERATURE AVAILABLE FROM SIGNET ARMORLITE:

- Swim / Ski / Dive Patient Pamphlet #480-554 (25/pk)
- Swim / Ski / Dive Sell Sheet (includes fitting suggestion/formulas for sports eyewear) #360-616
- Press-On Info / Instruction Booklet #358-925

3M Authorized Distributor
Medical Specialties Products

Please Recycle

3M and Press-on are trademarks of 3M.
©3M 2022. P#877488 Revised: July 2022


Signet Armorlite, Inc.
Phone: 800-759-0075
www.signetarmorlite.com