

Signia 1.50 Trifocal 8X35 76mm
Base Curve Selection Chart

SPHERE	CYLINDER																SPHERE	
	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75		-4.00
5.75																		5.75
5.50																		5.50
5.25																		5.25
5.00																		5.00
4.75																		4.75
4.50			8.75 Base															4.50
4.25																		4.25
4.00																		4.00
3.75																		3.75
3.50																		3.50
3.25																		3.25
3.00																		3.00
2.75																		2.75
2.50																		2.50
2.25																		2.25
2.00																		2.00
1.75																		1.75
1.50																		1.50
1.25																		1.25
1.00		6.25 Base																1.00
0.75																		0.75
0.50																		0.50
0.25																		0.25
0.00																		0.00
-0.25																		-0.25
-0.50																		-0.50
-0.75																		-0.75
-1.00																		-1.00
-1.25																		-1.25
-1.50																		-1.50
-1.75																		-1.75
-2.00																		-2.00
-2.25																		-2.25
-2.50			4.50 Base															-2.50
-2.75																		-2.75
-3.00																		-3.00
-3.25																		-3.25
-3.50																		-3.50
-3.75																		-3.75
	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00	